

# **Astral**<sup>™</sup>series

### Life support ventilators

## Astral with PS/SVt

Pressure Support with Safety Tidal Volume (PS/SVt) mode is a Volume Assured Pressure Support (VAPS) mode often used to treat patients with chronic obstructive pulmonary disease (COPD) and neuromuscular disease (NMD). The PS/SVt mode offers users a volume guarantee with every breath. This mode also provides patients with flexibility in using a traditional vented mask with standard tubing (6 ft). PS/SVt is now available on the Astral™ series leak circuit, which enhances leak compensation with Vsync technology for patients with persistent mask or oral leak.



### Conversion recommendations

Trilogy™ AVAPS AE settings	Astral PS/SVt conversion example	
Tidal Volume (Vte) Example: 500 mL	Safety Vt = Tidal Volume  Example: 500 mL	
Maximum Pressure Support (PS Max) $Example: PS max = 25 cmH_2O$	PS Max = PS Max  Example: Set at 25 cmH <sub>2</sub> O	
Minimum Pressure Support (PS Min)  Example: PS Min = 6 cmH <sub>2</sub> O	PS = PS Min  Example: Set at 6 cm $H_2O$	
Expiratory Positive Airway Pressure Maximum (EPAP Max)  Example: EPAP = 15 cmH <sub>2</sub> O	N/A	
Expiratory Positive Airway Pressure Minimum (EPAP Min) $Example: EPAP = 5 cmH_2O$	PEEP = EPAP  Example: Set PEEP at 5 cmH <sub>2</sub> O	
Respiratory Rate  Example: BPM = Auto*	Respiratory Rate = RR  Example: 10 bpm*	

<sup>\*</sup> There is no Auto backup rate in PS/SVt . Respiratory Rate should be set per physician's orders.

#### Pressure Support with Safety Tidal Volume (PS/SVt)

PS/SVt is a commonly used mode in COPD patients. PS/SVt is similar to Trilogy AVAPS, as AVAPS is also a VAPS mode.

#### **Primary PS/SVt settings**

Below are primary settings to be prescribed for PS with Safety Vt:

PS: Pressure Support. Set as the minimum pressure support requirement, also equivalent to Trilogy PS Min.

PS Max: Pressure Support Maximum. Set as the maximum pressure support allowable/limit, equivalent to Trilogy PS Max.

PEEP: Positive End-Expiratory Pressure, also equivalent to Trilogy EPAP.

RR: Respiratory Rate to be set as a traditional back-up rate, equivalent to Trilogy Rate.

Safety Vt: Volume guarantee, set at 6-8 mL/kg IBW or equivalent to Trilogy Tidal Volume on an AVAPS prescription.



### Astral specifications

totiai opeemeatiene	Astral 100	Astral 150
Indications for use	Life support for patients over 5 kg*	Life support for patients over 5 kg*
Circuit		
Single circuit with leak	•	•
Single circuit with valve	•	•
Double circuit		•
Therapy		
Valve circuit therapy: CPAP, (A)CV, P(A)CV, P-SIMV, V-SIMV, PS	•	•
Leak circuit therapy: CPAP, (S)T, P(A)C, iVAPS*, PS	•	•
Manual breath		•
Sigh breath (recruitment)		•
Apnea ventilation	Volume and pressure	Volume and pressure
Preset programs	2	4
Pressure range and settings		
Delivered pressure (cm H <sub>2</sub> 0)	2–50 (leak circuit) 0–50 (valve circuit)	2–50 (leak circuit) 0–50 (valve circuit)
Delivered tidal volume (cm H <sub>2</sub> O)	Adult: 100 to 2500 mL Pediatric: 50 to 300 mL	Adult: 100 to 2500 mL Pediatric: 50 to 300 mL
Respiratory rate	Adult: Off, 2 to 50 bpm Pediatric: Off, 5 to 80 bpm	Adult: Off, 2 to 50 bpm Pediatric: Off, 5 to 80 bpm
Data		
SpO <sub>2</sub> monitoring	Optional	Optional
FiO <sub>2</sub> monitoring	Optional	Standard
On screen	30 days of summary data Real-time pressure and flow (25 Hz)	30 days of summary data Real-time pressure and flow (25 H
ResScan™ (via download)	365 days of summary data 7 days of detailed data	365 days of summary data 7 days of detailed data
AirView <sup>™</sup> (RCM required)	365 days of summary data 7 days of detailed data	365 days of summary data 7 days of detailed data
Alarms	A wide range of settable alarms are available, including $\text{FiO}_2$ and $\text{SpO}_2$ when connected.	
Dimensions	-	
Weight	3.2 kg	3.2 kg
Dimensions (L x W x H)	285 mm x 215 mm x 93 mm	285 mm x 215 mm x 93 mm
Internal battery	Lithium-lon battery, 14.4 V, 6.6 Ah, 95 Wh. Up to 8 hours of backup power	Lithium-Ion battery, 14.4 V, 6.6 Ah, 95 Wh. Up to 8 hours of backup power
Sound power level	43 ± 3 dBA***	43 ± 3 dBA***
Sound pressure level	35 ± 3 dBA**	35 ± 3 dBA**
IEC 60601-1 CLASSIFICATIONS: Medical electrical equipment – general requirements	Class II (double insulation) Type BF Continuous operation	

## $^{\ast}$ iVAPS therapy mode is indicated for patients weighing 30 kg and above $^{\ast\ast}$ As measured according to ISO 80601-2-12:2011

for basic safety and essential performance.

### Product and accessory codes





Device trolley 279001



Astral mobility bag 27917



Astral SlimFit mobility bag 27949



Astral table stand 27915



External battery 27919



Astral Remote Alarm II (does not include cable) 27903



Single circuit smooth bore 15 mm  $\mathsf{QuickConnect}^{\scriptscriptstyle\mathsf{TM}}\ \mathsf{circuit}$ 21312



ResMed Connectivity Module (requires rental agreement) 27201

See the ventilation accessories catalog for the full suite of ResMed accessories.





















Suitable for use with oxygen

























<sup>1</sup> Köhnlein T et al. Non-invasive positive pressure ventilation for the treatment of severe stable chronic obstructive pulmonary disease: a prospective, multicentre, randomised, controlled clinical trial. Lancet Respir Med 2014;2:698–705. 2 Jaye J et al. Eur Respir J 2012

See ResMed.com for ResMed locations worldwide. AirView, Astral, QuickConnect and ResScan are trademarks and/or registered trademarks of the ResMed family of companies. Trilogy is a trademark of Koninklijke Philips N.V., registered in the U.S. and other countries. For patent information, see ResMed.com/ip. Specifications may change without notice. Product availability may vary across regions. © 2018 ResMed Ltd. 10110668/1 2018-04