


Features designed for patient usability



- 1 **Quick fitting headgear** with easy adjusting notches on the top strap. Set once and forget!
- 2 **Two-point support** has minimal contact on the face for all-round vision plus 'lock and key' for easy headgear detachment.
- 3 **SoftEdge™ fabric back strap** is soft and comfortable, with quick-setting fasteners.
- 4 **All-soft nasal cushion** with a modified dual-wall design and Spring-Air Technology™ is comfortably soft to wear.
- 5 **Ultra smooth ball-joint** with tube swivel provides 360° of unimpeded movement.
- 6 **Whisper-quiet vents** distribute air away from you and your bed partner gently and quietly—just 25 dBA.*
- 7 **Cushion tab hinges** adapt to improve flexibility between the cushion and the headgear.
- 8 **Spring-flex tubing** minimises tube drag on the mask to provide a dependable seal.

* Declared dual numbered noise emission values (with uncertainty 3 dBA) in accordance with ISO 4871.e1.

RESMED



To get more out of your therapy, pair your Swift FX Nano with ResMed's S9™ Series of sleep therapy devices. Visit www.resmed.com.au for more information.



ResMed Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia.

Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA.
[EC REP] ResMed (UK) Ltd 96 Jubilee Ave Milton Park Abingdon Oxfordshire OX14 4RW UK.
See www.resmed.com for other ResMed locations worldwide. Spring Air is a trademark of ResMed Ltd. S9, SoftEdge and Swift are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. © 2013 ResMed Ltd. Specifications may change without notice. 1016419/1 2013-06

RESMED

Small in size,
big on comfort



Swift™ FX Nano
NASAL MASK

Swift™ FX Nano for Her
NASAL MASK

RESMED



Open, streamlined design

The Swift FX Nano nasal mask combines the open, streamlined design of a nasal pillows mask that many users love, with the natural breathing comfort of a nasal mask. Its two-point headgear ensures minimal contact on the face and fits quickly with minimal adjustment making set-up easier. So now you can easily use your mask while you wear glasses, read or watch television.

Slip into something compact and comfortable

Swift FX Nano is flexible during movement, so you can sleep in a number of different positions including side sleeping, while maintaining a comfortable seal and enjoying a high level of comfort. With its slip-on headgear and all-soft cushion, Swift FX Nano is simple to use.

Comfort makes a world of difference

Swift FX Nano is designed to deliver an excellent user experience, without compromising on the big things that matter... like comfort, fit and ease of use. It comes with soft wraps for added comfort and to help reduce facial marking. The headgear is quick and easy to fit, and you only need to set it once. Plus with just three main parts, wearing, cleaning and assembling your mask has never been easier.



Swift FX Nano **for Her**

Swift FX Nano for Her offers women all the benefits of Swift FX Nano along with a specially designed cushion and smaller headgear to offer a fit range specifically for women. Whether you wear your hair short or long, under or over the pink backstrap, it complements your personal style.

