



ResMed

Astral™ Mobility Bag

ResMed's Astral Mobility Bag has been designed to protect the Astral ventilator while giving your patients greater freedom and mobility.



Made specifically for ResMed's Astral ventilator, the Astral Mobility Bag has been built for freedom and convenience.

Its compact and purpose-built design provides effective protection for the device, as well as a neat and intelligent solution for patients and their carers when mobility is key.



Technical specifications

Nominal dimensions	420 mm x 320 mm x 190 mm (L x W x H) (16.5" x 12.6" x 7.5")
Weight (empty)	Approx. 2.3 kg (5.07 lb)
Water ingress protection	Astral and external battery protected against dripping water when tilted up to 15 degrees from specified orientation (IPX2).
Alarm volume range	45–75 dBA (in five steps) (when used inside the Astral Mobility Bag)

More mobility and freedom

Broad, comfortable shoulder straps make the bag easy to carry either by hand or as a backpack, while additional straps allow for easy mounting onto wheelchairs.

Bags of space

Specific compartments within the bag have been thoughtfully designed to make room for all elements of the system—including the ventilator, external battery, and power supply unit.

Accessibility

The transparent panel on the protective cover makes the device's interface fully accessible while still in the bag. Therapy can easily be checked, and preset programs can be switched while on the move.

Product codes

Astral Mobility Bag	27917
---------------------	--------------

Discover more about Astral
at ResMed.com/Astral