Actual Size: 140mm (Width) x 175mm (Height)

Ruler is for scaling purposes.

How to use the fitting template:

- 1. Please ensure the template is printed to scale
- 2. Cut out along dotted line
- 3. Follow directions on template to determine appropriate cushion size




AirFit[™] N10

Nasal Mask for Her

Swift[™] FX Nano Swift[™] FX Nano for Her

AirFit and Swift are trademarks of ResMed Ltd.
Swift is registered in the U.S. Patent and Trademark Office.
© 2014 ResMed Ltd. 1017456/2 2014-08





Sample

size

check

10

20

30

40

50

60

70

80

90

100

110

120

130

140

150

160

170